

Standard Practice for Removal of Uranium or Plutonium, or both, for Impurity Assay in Uranium or Plutonium Materials¹

This standard is issued under the fixed designation C1647; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ε) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This practice covers instructions for using an extraction chromatography column method for the removal of plutonium or uranium, or both, from liquid or digested oxides or metals prior to impurity measurements. Quantification of impurities can be made by techniques such as inductively coupled plasma mass spectrometry (ICP-MS), inductively coupled plasma atomic emission spectrometry (ICP-AES) or atomic absorption spectrometry (AAS.)

1.2 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.

1.3 This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

2. Referenced Documents

2.1 ASTM Standards:²

- C753 Specification for Nuclear-Grade, Sinterable Uranium Dioxide Powder
- C757 Specification for Nuclear-Grade Plutonium Dioxide Powder, Sinterable
- C776 Specification for Sintered Uranium Dioxide Pellets
- C787 Specification for Uranium Hexafluoride for Enrichment
- C788 Specification for Nuclear-Grade Uranyl Nitrate Solution or Crystals

C859 Terminology Relating to Nuclear Materials

C996 Specification for Uranium Hexafluoride Enriched to Less Than 5 $\%\ ^{235}\text{U}$

- C1168 Practice for Preparation and Dissolution of Plutonium Materials for Analysis
- C1287 Test Method for Determination of Impurities in Nuclear Grade Uranium Compounds by Inductively Coupled Plasma Mass Spectrometry
- C1347 Practice for Preparation and Dissolution of Uranium Materials for Analysis
- C1432 Test Method for Determination of Impurities in Plutonium: Acid Dissolution, Ion Exchange Matrix Separation, and Inductively Coupled Plasma-Atomic Emission Spectroscopic (ICP/AES) Analysis
- C1517 Test Method for Determination of Metallic Impurities in Uranium Metal or Compounds by DC-Arc Emission Spectroscopy
- D1193 Specification for Reagent Water

3. Terminology

3.1 For definitions of terms used in this test method but not defined herein, refer to Terminology C859.

4. Summary of Practice

4.1 An aliquot of liquid sample or dissolved solid sample is adjusted as needed to 8M nitric acid for plutonium/uranium removal using extraction chromatography. Uranium and plutonium are retained on the resin and trace impurities are collected in the column effluent. The impurities can be measured by a variety of techniques.

5. Significance and Use

5.1 This practice can be used to separate uranium or plutonium, or both, prior to the impurity analysis by various techniques. The removal of uranium and plutonium prior to quantification can improve the detection limits by minimizing the signal suppression caused by uranium or plutonium when using ICP techniques. Detection limits of $\sim 1-10$ part-perbillion (PPB) may be obtainable by matrix removal. Also, removal of the uranium and plutonium may allow the impurities analysis to be performed on a non-glove box enclosed instrument.

5.2 Other test methods exist to determine impurities in uranium or plutonium. Test Method C1517 is able to determine many impurities in uranium at detection levels of \sim 1–10

¹ This practice is under the jurisdiction of ASTM Committee C26 on Nuclear Fuel Cycle and is the direct responsibility of Subcommittee C26.05 on Methods of Test.

Current edition approved Jan. 1, 2013. Published January 2013. Originally approved in 2006. Last previous edition approved in 2006 as C1647 – 06. DOI: 10.1520/C1647-13.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

part-per-million (ppm) by DC-Arc Spectrometry. Test Method C1287 is able to determine impurities in uranium at detection levels of ~100 ppb by ICP-MS. Test Method C1432 provides an alternative technique to remove plutonium by ion exchange prior to analysis of the impurities by ICP-AES.

5.3 This practice can be used to demonstrate compliance with nuclear fuel specifications, for example, Specifications C753, C757, C776, C787, C788, and C996.

6. Interferences

6.1 Zirconium, hafnium, niobium, and tantalum are retained on the diamyl, amylphosphonate resin unless hydrofluoric acid is added to the nitric acid load/rinse solution. The addition of hydrofluoric acid to recover zirconium, hafnium, niobium, and tantalum reduces uranium and plutonium retention. For this reason, hydrofluoric acid levels must be minimized (typically <0.05M HF) in the load/rinse solution to prevent uranium or plutonium, or both, from eluting from the column into the trace metal fraction.

7. Apparatus

7.1 *Large Columns*, >13 mL capacity (inner diameter = 1.5 cm has been found acceptable) and reservoirs.

7.2 Plastic Collection Tubes, 50 mL.

7.3 Column Rack, used for gravity flow systems.

7.4 Polyethylene Frits for columns, 20 µm.

7.5 *Vacuum Box*—The use of a vacuum-assisted flow system permits the use of higher eluent flow rates. Gravity flow systems may be used instead.

8. Reagents

8.1 *Purity of Reagents*—Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that all reagents shall conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society where such specifications are available.³ Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination. High purity acids may be used to reduce reagent blanks and to achieve lower detection limits.

8.2 *Purity of Water*—Unless otherwise indicated, references to water shall be understood to mean reagent water as defined in Specification D1193 as Type I.

8.3 *Hydrofluoric Acid* (28 M)—Concentrated hydrofluoric acid (sp gr 1.2).

8.4 Nitric Acid (16 M)—Concentrated HNO₃ (sp gr 1.42).

8.5 Nitric Acid Solution (8 M)—Add 500 mL of concentrated nitric acid (sp gr 1.42) to 300 mL of water and dilute to 1 L with water.

8.6 *Nitric Acid* (8M)—*Hydrofluoric Acid* (0.05M)—Add 500 mL concentrated nitric acid and 1.8 mL concentrated hydrofluoric acid to 250 mL water and dilute to 1 liter with water.

8.7 *Diamyl, Amylphosphonate Resin*,⁴ 50–100 μ m particle size resin for use with vacuum-assisted flow systems; 100 to 150 μ m if using gravity flow systems.

8.8 *Polymethacrylate Resin*, 5 100–150 μ m particle size resin.

9. Hazards

9.1 Refer to the laboratory's chemical hygiene plan and other applicable guidance for handling chemical and radioactive materials and for the management of radioactive, mixed, and hazardous waste.

9.2 Hydrofluoric acid is a highly corrosive acid that can severely burn skin, eyes, and mucous membranes. Hydrofluoric acid is similar to other acids in that the initial extent of a burn depends on the concentration, the temperature and the duration of contact with the acid. Hydrofluoric acid differs from other acids because the fluoride ion readily penetrates the skin causing destruction of deep tissue layers. Unlike other acids that are rapidly neutralized, hydrofluoric acid reactions with tissue may continue for days if left untreated. Due to the serious consequences of hydrofluoric acid burns, prevention of appropriate laboratory controls (hoods) and wearing adequate personal protective equipment to protect from skin and eye contact is essential.

10. Procedure

10.1 Column Preparation and Conditioning:

Note 1—If prepackaged columns are available from the resin supplier, skip to 10.1.3.

Note 2—Each column shall contain 10 mL diamyl, amylphosphonate resin. This resin amount can remove approximately 250 mg combined plutonium/uranium to enable impurity assay by ICP-AES, ICP-MS, or AAS.

10.1.1 Add approximately 3 mL of polymethacrylate resin to the bottom of a column. Rinse the column walls with water. Place a frit on top of the resin.

10.1.2 Add approximately 10 mL of diamyl, amylphosphonate resin to the column. Rinse the column walls with water and place another frit on top of the resin.

³ Reagent Chemicals, American Chemical Society Specifications, American Chemical Society, Washington, DC. For Suggestions on the testing of reagents not listed by the American Chemical Society, see Analar Standards for Laboratory Chemicals, BDH, Ltd., Poole, Dorset, U.K., and the United States Pharmacopeia and National Formulary, U.S. Pharmacopeial Convention, Inc. (USPC), Rockville, MD.

⁴ The only suppliers of the UTEVA Resin known to the committee at this time are Eichrom Technologies LLC, Lisle, IL, USA and Triskem International, Bruz, France. If you are aware of alternative suppliers, please provide this information to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee,¹ which you may attend. This resin is described in Horwitz, E. P., et al, "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography," Analytica Chimica Acta, 266, 1992, pp. 25-37.

⁵ The only suppliers of the Pre–filter material known to the committee at this time are Eichrom Technologies LLC, Lisle, IL, USA and Triskem International, Bruz, France. If you are aware of alternative suppliers, please provide this information to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee,¹ which you may attend.

10.1.3 Place the columns on a vacuum box if using a vacuum-assisted flow system or in a column rack if using a gravity flow system.

Note 3—If not using a vacuum box, the larger particle size resin should be used.

10.1.4 Add 10 mL 8M nitric acid to each column to condition the resin. Turn on vacuum and adjust flow rate to 2–4 mL/min if using a vacuum-assisted flow system. Allow each column to drain completely and turn off the vacuum if using the vacuum assisted flow system.

10.1.5 Discard the eluted conditioning solutions.

10.2 Column Loading and Separation:

10.2.1 Procedures for dissolution of plutonium and uranium materials are found in Practices C1168 and C1347.

10.2.2 Adjust the solution resulting from dissolution so that the nitric acid concentration is 8M.

10.2.3 Record the sample aliquot weight (g) or volume (mL) and sample pre-dilution factor (weight/volume or volume/volume).

10.2.4 Place clean, labeled 50 mL plastic collection tubes below each column in the vacuum extraction system.

10.2.5 Pipet 10 mL of sample aliquot prepared in 10.2.2 into the appropriate column. (the resin capacity is 250 mg U and Pu combined. If the concentration of the sample is greater than 25 mg/mL, then less than 10 mL should be added to avoid exceeding the resin capacity).

10.2.6 Turn on vacuum and adjust flow rate to 2–4 mL/min if using a vacuum-assisted flow system. Allow each column to drain completely.

10.2.7 Pipet 10 mL 8M nitric acid to each column and allow to drain completely.

10.2.8 Repeat step 10.2.7.

Note 4—This second column elution step may be performed with hydrofluoric added (20 mL 8M $\text{HNO}_3 - 0.05\text{M}$ HF) to completely elute zirconium, hafnium, tantalum and niobium (see 6.1). This is done in a second elution rinse to minimize hydrofluoric acid in the initial column elution so that silicon background at the ICP-AES is minimized and silicon can be measured in the first elution. A small amount of zirconium, hafnium, tantalum and niobium may be in the first elution so results from both elutions must be added to get a final result. If silicon is not of interest,

this step can be avoided by making the 8M nitric acid load solution approximately 0.05M hydrofluoric acid.

10.2.9 After columns have drained completely, increase the vacuum to remove residual liquid from columns, if using the vacuum-assisted flow system.

10.2.10 Turn off vacuum if using the vacuum-assisted flow system and remove collection tubes from vacuum system. Adjust collection tube volume to 30 mL with 8M nitric acid, if needed.

10.2.11 Submit collected fractions for analysis by the selected technique such as ICP-MS, ICP-AES, or AAS.

10.3 While not strictly within the scope of this practice, plutonium or uranium, or both, retained on the column may be eluted at this point, if desirable for the purposes of the laboratory performing this analysis. Guidance for accomplishing this may be found in Horwitz, et al.⁴

11. Calculations

11.1 Calculate the Metal Impurity, *I*:

$$I = (S - B) \times D \tag{1}$$

where:

S = metallic impurity measured in final solution, $\mu g/mL$,

- B = blank measurement for metallic impurity in final solution, μ g/mL, and
- D = sample dilution factor.

$$D = \frac{SV \times P}{SA} \tag{2}$$

where:

- SV = column removal volume, mL (typically 30 mL, step 10.2.10),
- SA = sample aliquot used in column separation, mL (typically 10 mL, step 10.2.5), and

P = pre-dilution factor (step 10.2.3).

12. Keywords

12.1 AAS; extraction chromatography; ICP-AES; ICP-MS; impurities; plutonium; resin; separation; uranium

ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org). Permission rights to photocopy the standard may also be secured from the ASTM website (www.astm.org/COPYRIGHT/).