

Standard Test Method for Analysis of High Manganese Steel by Spark Atomic Emission Spectrometry¹

This standard is issued under the fixed designation E2209; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This test method provides for the analysis of high manganese steel by spark atomic emission spectrometry for the following elements in the ranges shown:

| Elements Conposition Ra | |
|-------------------------|---------------|
| Aluminum (AI) | 0.02 to 0.15 |
| Carbon (C) | 0.3 to 1.4 |
| Chromium (Cr) | 0.25 to 2.00 |
| Manganese (Mn) | 8.0 to 16.2 |
| Molybdenum (Mo) | 0.03 to 2.0 |
| Nickel (Ni) | 0.05 to 4.0 |
| Phosphorus (P) | 0.025 to 0.06 |
| Silicon (Si) | 0.25 to 1.5 |
| | |

Note 1—The ranges represent the actual levels at which this method was tested.² These composition ranges can be extended to higher compositions by the use of suitable reference materials. Sulfur is not included because differences in results between laboratories exceeded acceptable limits at all analyte levels.

1.2 This test method may involve hazardous materials, operations, and equipment. This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

2. Referenced Documents

2.1 ASTM Standards:³

A128/A128M Specification for Steel Castings, Austenitic Manganese

E135 Terminology Relating to Analytical Chemistry for Metals, Ores, and Related Materials

E305 Practice for Establishing and Controlling Atomic

Emission Spectrochemical Analytical Curves

E353 Test Methods for Chemical Analysis of Stainless, Heat-Resisting, Maraging, and Other Similar Chromium-Nickel-Iron Alloys

E406 Practice for Using Controlled Atmospheres in Spectrochemical Analysis

E1019 Test Methods for Determination of Carbon, Sulfur, Nitrogen, and Oxygen in Steel, Iron, Nickel, and Cobalt Alloys by Various Combustion and Fusion Techniques

E1059 Practice for Designating Shapes and Sizes of Nongraphite Counter Electrodes (Withdrawn 2013)⁴

E1329 Practice for Verification and Use of Control Charts in Spectrochemical Analysis

E1601 Practice for Conducting an Interlaboratory Study to Evaluate the Performance of an Analytical Method

E1806 Practice for Sampling Steel and Iron for Determination of Chemical Composition

2.2 Other Document:

ASTM MNL 7 ASTM Manual on Presentation of Data and Control Chart Analysis, 8th Edition, 2010.

3. Terminology

3.1 For definition of terms used in this method, refer to Terminology E135.

4. Summary of Test Method

4.1 A controlled discharge is produced between the flat surface of the specimen and the counter electrode. The radiant energies of selected analytical lines are converted into electrical energies by photomultiplier tubes and stored on capacitors. This discharge is terminated after a fixed integration time. At the end of the exposure period, the charge on each capacitor is measured and converted to mass fraction percent.

5. Significance and Use

5.1 The chemical composition of high manganese steel alloys must be determined accurately to ensure the desired metallurgical properties. This procedure is suitable for manufacturing control and inspection testing.

¹ This test method is under the jurisdiction of ASTM Committee E01 on Analytical Chemistry for Metals, Ores, and Related Materials and is the direct responsibility of Subcommittee E01.01 on Iron, Steel, and Ferroalloys.

Current edition approved Nov. 15, 2013. Published January 2013. Originally approved in 2002. Last previous edition approved in 2009 as E2209 – 02 (Reapproved 2006) $^{\rm e2}$. DOI: 10.1520/E2209-13.

² Supporting data have been filed at ASTM International Headquarters and may be obtained by requesting Research Report RR:E01-1035.

³ For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

⁴ The last approved version of this historical standard is referenced on www.astm.org.

6. Interferences

6.1 Interferences may vary with spectrometer design and excitation characteristics. Direct spectral interferences may be present on one or more of the wavelengths listed in a method. Frequently, these interferences may be determined and proper corrections made by the use of various reference materials. The composition of the sample being analyzed should match closely the composition of one or more of the reference materials used to prepare and control the calibration curve that is employed. Alternatively, mathematical corrections may be used to solve for interelement effects. Various mathematical correction procedures are commonly utilized. Any of these are acceptable that will achieve analytical accuracy equivalent to that provided by this method.

7. Apparatus

- 7.1 Sample Preparation Equipment:
- 7.1.1 Sample Mold, to produce chilled cast samples approximately 38 mm (1½ in.) in diameter that are homogeneous, free of voids or porosity in the region to be excited, and representative of the material to be analyzed. Refer to Practice E1806 for steel sampling procedures.
- 7.1.2 *Immersion Sampler*, to take a sample from the bath or from the metal stream when pouring can be used. The sample should produce a sample of the same dimensions as listed in 7.1.1.
- 7.1.3 Surface Grinder or Sander With Abrasive Belts or Disk, capable of providing a flat uniform surface on the reference materials and specimens. The following table shows the various methods of sample preparation used in the Inter-Laboratory Study (ILS):

Type of Grinding Preparation Grinding Medium

Belt and/or Disk Aluminum Oxide, Zirconium Oxide 36 to 180

Grit of Grinding Medium

Note 2—Silicon carbide grinding medium may be used but it was not utilized by the laboratories in the Inter-Laboratory Study (ILS).

- 7.2 Excitation Source, capable of providing a triggered capacitor discharge having the source parameters meeting the requirements of 11.1.
- 7.3 Excitation Stand, suitable for mounting in optical emission alignment, a flat surface for the specimen in opposition to a counter electrode. This stand shall provide an atmosphere of argon. The electrode and argon are described in 8.1 and 8.2.
- 7.4 Spectrometer, having sufficient resolving power and linear dispersion to separate clearly the analytical lines from other lines in the spectrum of a specimen in the spectral region 170.0 nm to 450 nm. The spectrometer shall have a dispersion of at least 2 nm/mm and a focal length of at least 0.5 m. Gas purged spectrometers are an alternative to vacuum systems.
- 7.5 Measuring System, consisting of photomultiplier tubes having individual voltage adjustment, capacitors on which the output of each photomultiplier tube is stored and an electronic system to measure voltages on the capacitors either directly or indirectly, and the necessary switching arrangements to provide the desired sequence of operation.

- 7.6 *Vacuum Pump*, if required, capable of maintaining a vacuum of approximately 3 Pa. There are some equipment manufactures that will purge the optical portion of the spectrometer with argon or other inert gas rather than pull a vacuum on the optics. Either vacuum optics or purged optics are required to determine carbon and phosphorus in this method.
- 7.7 Gas System, consisting of an argon supply, a pressure regulator, and a gas flow meter. Automatic sequencing shall be provided to actuate the flow of argon at a given flow rate for a given time interval and to start the excitation at the end of the required flush period. The gas system shall be in accordance with Practice E406.

8. Reagents and Materials

- 8.1 *Argon*, either gaseous or liquid, must be of sufficient purity to permit proper excitation of the analytical lines of interest. Argon of 99.998 % purity has been found satisfactory. Refer to Practice E406.
- 8.2 Counter Electrode—A Tungsten or Thoriated Tungsten rod ground to a 15°, 30°, 45° or 90° angle conical tip, which conforms to Practice E1059, was found satisfactory. Other material may be used provided it can be shown experimentally that equivalent precision and bias are obtained.

9. Reference Materials

- 9.1 *Certified Reference Materials*, for high manganese steel are commercially available.
- 9.2 Calibrants shall be certified reference materials from recognized certification agencies. They shall cover the composition ranges of the elements to be determined and shall include

TABLE 1 Wavelengths

| Element | Wavelength | Line | Possible | | |
|--------------------------|------------|----------------|----------------------------|--|--|
| Element | (nm) | Classification | Interferences ^A | | |
| Aluminum | 394.4 | I | V, Mn, Mo | | |
| | 396.152 | 1 | Мо | | |
| Carbon | 193.09 | 1 | Al | | |
| Chromium | 298.92 | II | Mn, V, Ni, Nb, Mo | | |
| | 267.72 | II | Mn, Mo, V | | |
| | 425.435 | I | | | |
| Iron (Internal Standard) | 273.07 | I | | | |
| | 271.44 | II | | | |
| Manganese | 263.81 | II | | | |
| | 290.02 | II | | | |
| | 293.31 | II | Cr | | |
| Molybdenum | 202.03 | II | | | |
| | 263.876 | II | | | |
| | 281.61 | II | Al, Mn | | |
| | 386.41 | I | V, Cr | | |
| Nickel | 231.60 | II | Co, Ti | | |
| | 218.54 | II | | | |
| | 352.45 | I | | | |
| | 341.476 | 1 | | | |
| Phosphorus | 178.29 | 1 | Мо | | |
| Silicon | 212.41 | I | | | |
| | 288.16 | 1 | Mo, Cr, W | | |
| | 251.61 | I | Fe, V | | |
| Sulfur | 180.73 | I | Mn | | |

^A Interferences are dependent upon instrument design, and excitation conditions, and those listed require confirmation based upon specimens designed to demonstrate interferences. This standard method does not purport to address the interferences that these lines may have. Care should be taken to address the interferences when calibrating the instrument.



all of the specific types of alloys being analyzed. The calibrants shall be homogeneous and free of voids and porosity. The metallurgical history of the calibrants should be similar to that of the specimens being analyzed. Refer to Test Methods E353 and E1019 for chemical analysis of high manganese steel alloys.

9.2.1 In selecting calibrants, use caution with compositions that are unusual. One element may influence the radiant energy of another element. Tests should be made to determine if interrelations exist between elements in the calibrants.

10. Preparation of Calibrants and Specimens

10.1 Rough grind, either wet or dry, with a coarse grinding belt or disk. The final grind of the specimen must be the same grit as the calibrants. Dry the specimens, if wet, for proper excitation in the argon atmosphere. Ensure that the specimens are homogeneous and free from voids and pits in the region to be excited. Cast specimens from molten metal into a suitable mold and cool. Immersion and stream samplers are also suitable for use. Prepare the surface of the specimens and reference materials in a similar manner.

11. Excitation and Instrument Parameters

11.1 *Electrical Parameters*—Electrical parameters within the following ranges were found acceptable.

11.2 Spectrometer Configurations:

Spectrometer Parameters

 Focal Length
 0.5 m to 1.2 m

 Dispersion
 0.5 nm/mm to 2.16 nm/mm

 Vacuum
 1 Pa to 25 Pa

11.3 Excitation Conditions:

 Flush Time
 2 s to 5 s

 Preburn
 10 s to 30 s

 Exposure
 5 s to 20 s

11.4 *Initiation Circuit*—The initiator circuit parameters shall be adequate to uniformly trigger the capacitor discharge. The values for these parameters will vary with the instrument. Normal values found to be adequate are listed as follows:

Capacitance (d-c charged) :F 1.2 Inductance, :H residual Resistance, Σ residual Potential, V 425

- 11.4.1 *Other Electrical Parameters*—Excitation units, on which the precise parameters given in 11.1 and 11.4 are not available, may be used provided that it can be shown experimentally that equivalent precision and accuracy are obtained.
- 11.5 *Electrode System*—Insert the counter electrode in the lower electrode position. Adjust the analytical gap to 3 mm, 4 mm, 5 mm, or 7 mm depending on the manufacturer's recommendations for that particular instrument.

11.6 *Discharge Source*—Most capacitor discharge sources in today's spectrometers are either the directional self-initiating capacitor discharge type or a triggered capacitor discharge type.

12. Preparation of Instrumentation

12.1 Prepare the spectrometer in accordance with the manufacturer's instructions.

Note 3—It is not within the scope of this method to prescribe all details of equipment to be used. Equipment varies between laboratories.

13. Calibration, Standardization, and Verification

- 13.1 Calibration—Using the conditions given in 11.3, excite the calibrants and potential drift correction samples in a random sequence, bracketing these burns with excitations of any materials intended for use as verifiers. (A verifier may be used as a calibrant even though it is burned only as a verifier.) There should be at least five calibrants for each element, spanning the required composition range. Make replicate excitations in accordance with 14.2. Using the averages of the data for each point, determine analytical curves as described in Practice E305 or the instrument manufacturer's instructions..
- 13.2 Standardization—Following the manufacturer's recommendations, standardize on an initial setup and anytime that is known or suspected that readings have shifted. Make the necessary corrections by applying mathematical corrections. Standardization shall be done anytime verifications indicate that readings have gone out of statistical control.
- 13.3 *Verification*—Shall be done at least at the beginning of any analytical work. Analyze verifiers in replicate to confirm that they read within expected confidence interval, as defined in 13.4. The replication shall be the same as directed in 14.2.
- 13.3.1 Check the verification after standardizing. If confirmation is not obtained, standardize again and/or investigate why confirmation is not obtained. Standardization is confirmed if the results are within two standard deviations from the mean of the standard.
- 13.3.2 Repeat the verification at least every 4 h or if the instrument has been idle for more than 1 h. If readings are not in conformance, repeat the standardization.
- 13.4 The confidence interval will be established from observations of the repeatability of the verifiers and determining the confidence interval for some acceptable confidence level or by establishing the upper and lower limit of a control chart in accordance with PracticeE1329 or ASTM MNL 7.

14. Measurement Procedure for Excitation and Radiation

- 14.1 Use the electrical conditions given in Section 11.
- 14.2 Place the freshly-surfaced specimen on the excitation stand. Position it to effect a gas-tight seal and adequate argon flushing. Position the specimen so that there will be a uniform pattern of burns around its face. For example, a disk-shaped specimen should have a ring of burn marks around its outer edge and approximately 6 mm (0.25 in.) from the edge. Avoid burning the center of cast specimens where quench cracks and segregation are more likely to be present. Ensure there is a good electrical ground between the specimen and the specimen

holder device. Cool the specimen after two burns to prevent overheating, if required. One way of cooling the sample without water is to place it on a solid piece of copper to act as a heat sink. Examine the specimen after each burn to evaluate the quality of the excitation. Cracks, voids, pits, moisture, or inclusions will limit the sampling and accuracy of a determination. Successive burns shall be sufficiently separated so that the burn patterns do not overlap.

14.3 Replicate Excitations—As a minimum, make duplicate burns on each specimen. Average the duplicate readings for each element if their difference does not exceed twice the established standard deviation for that element. Most instruments today will print out a statistical analysis of the burns for the analysis. If their difference exceeds this value, analyze the specimens two more times. Examine the data to determine if there is an outlier. If there is an outlier, average the other three results. If there are two outliers, average all four burns or repeat the analysis. (Two outliers may indicate sample inhomogeneity).

15. Calculation of Results

15.1 Average the readings obtained for each specimen. If the readout is not in direct mass fraction units, use this intensity to obtain the mass fractions from the curves.

16. Precision and Bias⁵

16.1 *Precision*—The precision of this test method was determined by submitting 7 unknown samples to 9 different laboratories. The interlaboratory testing was conducted in accordance with Practice E1601. Using reference materials supplied with the unknowns, calibration curves were established on the equipment and the corresponding 7 unknown samples were analyzed against the established curves. Precision data are found in Table 2.

16.2 *Bias*—The difference between the certified composition and the average composition determined by this test method indicates the bias encountered in the interlaboratory test. Bias data is found in Table 2.

17. Keywords

17.1 spark atomic emission spectrometry; high manganese steel

⁵ Supporting data have been filed at ASTM International Headquarters and may be obtained by requesting Research Report RR:E01-1031.



TABLE 2 Precision and Bias Data

| Certified Composition, % | Average Composition, % | No. of Labs | Min SD (s _M , E1601) | Reproducibility SD (s _R , E1601) | Reproducibility Index (R, E1601) | R_{rel} % | Material Test Number |
|--------------------------|------------------------|-------------|------------------------------------|--|----------------------------------|-------------|-------------------------|
| Aluminum | | | | | | | |
| 0.01 | 0.0044 | 5 | 0.0005 | 0.0037 | 0.0103 | 237 | 7 |
| 0.004 | 0.0047 | 7 | 0.0005 | 0.0034 | 0.0096 | 203 | 5 |
| 0.042 | 0.0470 | 6 | 0.0012 | 0.0024 | 0.0067 | 14.2 | 3 |
| 0.103 | 0.1073 | 8 | 0.0012 | 0.0056 | 0.0158 | 14.7 | 6 |
| 0.14 | 0.1220 | 6 | 0.0018 | 0.0131 | 0.0367 | 30.0 | 4 |
| Carbon | | | | | | | |
| 0.303 | 0.338 | 9 | 0.00425 | 0.1070 | 0.2997 | 88.0 | 1 |
| 0.79 | 0.804 | 9 | 0.00562 | 0.0258 | 0.0722 | 8.99 | 7 |
| 0.82 | 0.813 | 9 | 0.00802 | 0.0132 | 0.0371 | 4.56 | 6 |
| 0.92 | 0.908 | 7 | 0.00707 | 0.0121 | 0.0338 | 3.72 | 3 |
| | | 9 | | | | | |
| 1.24 | 1.241 | | 0.00934 | 0.0165 | 0.0462 | 3.72 | 5 |
| 1.22 | 1.263 | 9 | 0.00804 | 0.0244 | 0.0682 | 5.40 | 2 |
| 1.4 | 1.395 | 7 | 0.01028 | 0.0163 | 0.0457 | 3.28 | 4 |
| Chromium | | | | | | | |
| 0.27 | 0.2770 | 9 | 0.0036 | 0.0214 | 0.0600 | 21.7 | 7 |
| 0.474 | 0.5098 | 9 | 0.0064 | 0.0568 | 0.1589 | 31.2 | 1 |
| 0.56 | 0.5439 | 9 | 0.0037 | 0.0234 | 0.0655 | 12.0 | 5 |
| 0.62 | 0.5842 | 9 | 0.0070 | 0.0286 | 0.0799 | 13.7 | 2 |
| 1.02 | 0.9621 | 7 | 0.0094 | 0.0741 | 0.2075 | 21.6 | 4 |
| 1.45 | 1.446 | 7 | 0.0094 | 0.0539 | 0.1509 | 10.4 | 3 |
| 1.93 | 1.924 | 9 | | 0.0296 | 0.1309 | 4.30 | 6 |
| | 1.924 | 9 | 0.0165 | 0.0290 | 0.0828 | 4.30 | Ö |
| Manganese | | _ | | | | | _ |
| 8.28 | 8.115 | 7 | 0.0886 | 0.1896 | 0.5310 | 6.54 | 4 |
| 9.25 | 8.990 | 9 | 0.0749 | 0.3416 | 0.9566 | 10.6 | 7 |
| 12.4 | 12.409 | 8 | 0.0562 | 0.2527 | 0.7076 | 5.70 | 2 |
| 13.6 | 13.582 | 8 | 0.0819 | 0.1375 | 0.3849 | 2.83 | 6 |
| 13.5 | 13.584 | 9 | 0.0927 | 0.1426 | 0.3992 | 2.94 | 5 |
| 15.6 | 15.784 | 9 | 0.1045 | 0.3754 | 1.051 | 6.66 | 1 |
| 16.1 | 16.152 | 7 | 0.1529 | 0.1819 | 0.5094 | 3.15 | 3 |
| | 10.132 | , | 0.1523 | 0.1019 | 0.3034 | 5.15 | 3 |
| Molybdenum | 0.0100 | 0 | 0.0010 | 0.0050 | 0.0100 | 00.0 | 0 |
| 0.023 | 0.0166 | 6 | 0.0019 | 0.0058 | 0.0163 | 98.2 | 2 |
| 0.01 | 0.0173 | 5 | 0.0016 | 0.0106 | 0.0298 | 172 | 4 |
| 0.035 | 0.0346 | 9 | 0.0111 | 0.0158 | 0.0444 | 128 | 6 |
| 0.078 | 0.0690 | 9 | 0.0130 | 0.0192 | 0.0538 | 77.9 | 5 |
| 0.6 | 0.5807 | 9 | 0.0180 | 0.0373 | 0.1044 | 18 | 5 3 |
| 1.55 | 1.543 | 7 | 0.0221 | 0.0849 | 0.2379 | 15.4 | 7 |
| 1.91 | 1.992 | 9 | 0.0367 | 0.2444 | 0.6843 | 34.3 | 1 |
| Nickel | | · · | 0.0007 | 0.2 | 0.00.0 | 00 | • |
| 0.05 | 0.0506 | 6 | 0.0028 | 0.0035 | 0.0097 | 19.2 | 3 |
| | | 8 | | | | | 2 |
| 0.118 | 0.1132 | | 0.0027 | 0.0145 | 0.0405 | 35.8 | 2 |
| 0.58 | 0.6193 | 8 | 0.0044 | 0.0185 | 0.0517 | 8.36 | 7 |
| 0.69 | 0.6923 | 8 | 0.0076 | 0.0130 | 0.0363 | 5.24 | 5 |
| 1.05 | 1.0555 | 8 | 0.0190 | 0.0375 | 0.1051 | 9.96 | 6 |
| 3.13 | 3.030 | 8 | 0.0302 | 0.0809 | 0.2266 | 7.48 | 1 |
| 3.9 | 3.966 | 6 | 0.0167 | 0.0740 | 0.2073 | 5.23 | 4 |
| Phosphorus | | | | | | | |
| 0.024 | 0.0242 | 9 | 0.0006 | 0.0028 | 0.0078 | 32.2 | 1 |
| 0.026 | 0.0247 | 7 | 0.0008 | 0.0016 | 0.0045 | 18.1 | 3 |
| | | _ | | | | | _ |
| 0.03 | 0.0299 | 9 | 0.0007 | 0.0023 0.0015 | 0.0063 | 21.2 | 7 |
| 0.036 | 0.0341 | 9 | 0.0010 | | 0.0043 | 12.6 | 6 |
| 0.033 | 0.0377 | 9 | 0.0006 | 0.0016 | 0.0046 | 12.2 | 2 |
| 0.04 | 0.0378 | 9 | 0.0010 | 0.0017 | 0.0048 | 12.6 | 5 |
| 0.06 | 0.0523 | 7 | 0.0014 | 0.0030 | 0.0085 | 16.2 | 4 |
| Silicon | | | | | | | |
| 0.25 | 0.2593 | 9 | 0.0028 | 0.0059 | 0.0164 | 6.34 | 5 |
| 0.374 | 0.3801 | 9 | 0.0040 | 0.0142 | 0.0397 | 10.4 | 2 |
| 0.46 | 0.4599 | 9 | 0.0031 | 0.0076 | 0.0214 | 4.65 | 6 |
| 0.67 | 0.6669 | 9 | 0.0046 | 0.0368 | 0.1030 | 15.4 | 1 |
| | 0.9164 | | | | | | |
| 0.9 | | 7 | 0.0078 | 0.0215 | 0.0602 | 6.57 | 3 |
| 0.98 | 0.9290 | 7 | 0.0098 | 0.0484 | 0.1355 | 14.6 | 4 |
| 1.65 | 1.477 | 9 | 0.0140 | 0.0475 | 0.1330 | 9.00 | 7 |
| Sulfur | | | | | | | |
| 0.001 | 0.0019 | 4 | 0.0006 | 0.0019 | 0.0052 | 272 | 2 |
| 0.002 | 0.0030 | 7 | 0.0005 | 0.0008 | 0.0022 | 73.3 | 1 |
| 0.005 | 0.0039 | 5 | 0.0007 | 0.0007 | 0.0020 | 51.9 | 5 |
| 0.003 | 0.0064 | 8 | 0.0007 | 0.0024 | 0.0068 | 107 | 7 |
| | | | | | | | |
| 0.012 | 0.0087 | 7 | 0.0007 | 0.0028 | 0.0080 | 91.4 | 3 |
| 0.014 | 0.0129 | 9 | 0.0011 | 0.0026 | 0.0074 | 57.0 | 6 |
| 0.03 | 0.0202 | 7 | 0.0011 | 0.0041 | 0.0115 | 56.8 | 4 |



ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org). Permission rights to photocopy the standard may also be secured from the ASTM website (www.astm.org/COPYRIGHT/).