

# Standard Test Method for Determining the Bacteria-Eliminating Effectiveness of Hygienic Handwash and Handrub Agents Using the Fingerpads of Adults<sup>1</sup>

This standard is issued under the fixed designation E2276; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon  $(\varepsilon)$  indicates an editorial change since the last revision or reapproval.

### INTRODUCTION

Hands can spread many types of pathogens directly (1)<sup>2</sup> or by transfer of such organisms to other surfaces and objects during casual contact (2,3). Therefore, regular and proper decontamination of hands by caregivers and food-handlers in particular is crucial for infection control. Hygienic hand antisepsis is meant to reduce the load of transient microflora on hands, thereby reducing the risk of disease transmission. Such reduction in the bacterial load may be due to a combination of bacterial inactivation and removal of viable bacteria from the skin. In this method the test bacterial suspension is placed on the thumb- and fingerpads of adults to simulate the contamination of hands with transient microflora, the inoculum on the fingerpads is allowed to dry and is then treated with test and control solutions. Since in each test all ten digits on any given subject can be used, the protocol permits the inclusion of the required controls and several replicates of the test formulation in the same trial.

### 1. Scope

- 1.1 This test method is designed to determine the activity of hygienic handwash and handrub (4) agents against transient bacterial flora on hands and is not meant for use with surgical hand scrubs or preoperative skin preps.
- 1.2 Performance of this procedure requires the knowledge of regulations pertaining to the protection of human subjects.<sup>3</sup>
- 1.3 The test method should be performed by persons with training in microbiology in facilities designed and equipped for work with infectious agents at biosafety level 2 (5).
- 1.4 The values stated in SI units are to be regarded as standard. No other units of measurement are included in this standard.
- 1.5 This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.
- <sup>1</sup> This test method is under the jurisdiction of ASTM Committee E35 on Pesticides, Antimicrobials, and Alternative Control Agents and is the direct responsibility of Subcommittee E35.15 on Antimicrobial Agents.
- Current edition approved April 1, 2010. Published May 2010. Originally approved in 2003. Last previous edition approved in 2003 as E2276 03  $^{\rm c1}$ . DOI: 10.1520/E2276-10.
- <sup>2</sup> The boldface numbers in parentheses refer to the list of references at the end of this standard.
  - <sup>3</sup> Federal Register, Vol 46, No. 17, Jan. 27, 1991.

#### 2. Referenced Documents

2.1 ASTM Standards:<sup>4</sup>

D1129 Terminology Relating to Water

E1115 Test Method for Evaluation of Surgical Hand Scrub Formulations

E1173 Test Method for Evaluation of Preoperative, Precatheterization, or Preinjection Skin Preparations

E1174 Test Method for Evaluation of the Effectiveness of Health Care Personnel Handwash Formulations

E1838 Test Method for Determining the Virus-Eliminating Effectiveness of Hygienic Handwash and Handrub Agents Using the Fingerpads of Adults

E2613 Test Method for Determining Fungus-Eliminating Effectiveness of Hygienic Handwash and Handrub Agents Using Fingerpads of Adults

### 3. Terminology

- 3.1 *Definitions*—For definitions of general terms used in this test method, refer to Terminology D1129.
  - 3.2 Definitions of Terms Specific to This Standard:
- 3.2.1 *active ingredient, n*—a substance added to a formulation specifically for the inhibition or inactivation of microorganisms.

<sup>&</sup>lt;sup>4</sup> For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

- 3.2.2 *dry control*, *n*—a control to determine the number of colony forming units of the test bacterium remaining viable after the initial drying of the inoculum on the skin.
- 3.2.3 handrub, n—a liquid or gel which is applied by rubbing to decontaminate lightly soiled hands between handwashings and generally do not require a post-treatment water rinse; such agents usually contain alcohol alone or with other active ingredients.
- 3.2.4 *hard water*; *n*—water with a standard hardness of calcium carbonate.
- 3.2.5 hygienic handwash agent, n—an agent generally used for handwashing by personnel in hospitals, other health-care facilities, day-care centers, nursing homes, and food-handling establishments to eliminate transient microorganisms from intact skin.
- 3.2.6 *input control*, *n*—a control to determine the number of colony forming units of the test bacterium placed on each digit.
- 3.2.7 *neutralization, n*—a process which results in quenching the antimicrobial activity of a test substance. This may be achieved through dilution of the test substance to reduce the antimicrobial activity, or through the use of chemical agents, called neutralizers, to eliminate antimicrobial activity.
- 3.2.8 *nonmedicated soap*, *n*—a soap or detergent that is mild to the skin and does not contain any antimicrobial chemicals.
- 3.2.9 *soil load*, *n*—a solution of one or more organic and/or inorganic substances added to the suspension of the test organism to simulate the presence of body secretions, excretions or other extraneous substances.
- 3.2.10 *test substance or test formulation*, *n*—a formulation which incorporates antimicrobial ingredients.
- 3.2.11 *test organism*, *n*—an applied inoculum of an organism that has characteristics which allow it to be readily identified. The test organism is used to simulate a transient topical microbial contaminant. It may also be referred to as a marker organism, bacterial simulant/surrogate or bacterial contaminant
- 3.2.12 *test vehicle*, *n*—the test formulation without an active ingredient.
- 3.2.13 *transient microbiota*, *n*—microorganisms from the environment that contaminate but do not normally colonize the skin.

### 4. Summary of Test Method

4.1 This test method is conducted on a group of adults who have provided informed consent and the skin of whose hands has been determined to be free from any apparent damage. Subjects are to refrain from using any products containing antimicrobial agents for one week prior to the test. A known volume of the test bacterial suspension is placed on a demarcated area on each fingerpad and the inoculum allowed to dry. The contaminated area then is exposed to the control (standard hard water) or test substance for the desired contact time and organisms remaining on the fingerpad are eluted and the eluates assayed for viable bacteria. Percent or log<sub>10</sub> reductions in the numbers of viable bacteria after treatment with the control and test substance are then determined. The fingerpad

method gives results that are comparable to those obtained using a whole-hand procedure (6). If two different formulations are being compared in the same test, one of them may be designated as a reference and used in place of the hard water control. If desired, one also may use tap water in parallel with the hard water control to determine the influence of water hardness on the test product's bacteria-eliminating activity.

# 5. Significance and Use

- 5.1 This *in vivo* procedure is designed to test the ability of hygienic handwash or handrub agents to eliminate selected types of bacteria from experimentally contaminated skin of the hands of adult subjects. Since the two thumbpads and all eight fingerpads can be used in any given test, it allows for the incorporation of an input control (two), control for viable bacteria remaining after the inoculum has been allowed to dry (two), bacteria eliminated after treatment with a control or reference solution (two), and up to four replicates to assess the bacteria-eliminating efficiency of the product under test. No more than  $100~\mu L$  of the test bacterial suspension is required to complete one test. The results of testing with this test method may form the basis for confirmatory tests using a suitable whole-hand test protocol, such as Test Method E1174.
- 5.2 Whereas this test method relates to testing with bacteria, it can be readily adapted to work with protozoa and bacteriophages. Similar methods for work with fungi (Test Method E2613) and viruses of human origin (Test Method E1838) are already ASTM standards.
- 5.3 Potentially infectious microorganisms left on hands after washing can be reduced further by drying the washed hands with paper, cloth, or warm air (7). A step for the drying of fingerpads after exposure to the control or test solution, therefore, has not been included to avoid bacterial removal by the drying process itself.
- 5.4 This test method is not meant for use with surgical hand scrubs (Test Method E1115) or preoperative skin preps (Test Method E1173).
- 5.5 The level of contamination with viable bacteria on each fingerpad after the drying of the inoculum should be five- to ten-fold higher than the product performance criterion required. For example, the titer in the dried inoculum on each fingerpad should be about  $10^5$  colony forming units of the test bacterium when a  $>10^4$  reduction is required under the conditions of this test method.

# 6. Equipment and Apparatus

- 6.1 *Colony Counter*—Any of several types may be used, for example, Quebec Colony Counter.
- 6.2 Freezers—A freezer at  $-20 \pm 2^{\circ}$ C is required for the storage of culture media. A second freezer at  $-70^{\circ}$ C or lower is required to store bacterial stocks.
- 6.3 *Handwashing Sink*—A sink of sufficient size to permit subjects to wash hands without touching hands to sink surface.
- 6.4 *Incubator*—Any incubator capable of maintaining a temperature of  $36 \pm 1^{\circ}$ C. *Serratia marcescens* requires incubation at  $25 \pm 2^{\circ}$ C for pigment formation.

- 6.5 Laminar Flow Cabinet—A Class II biological safety cabinet is required for this work. The procedures for the proper maintenance and use of such cabinets are given in Ref (2).
- 6.6 *Magnetic Stirrer and Magnets*—Large enough to hold a 5-L beaker or Erlenmeyer flask for preparing culture media or other solutions.
- 6.7 Membrane Filtration System—A membrane filtration system and membranes with a pore diameter of 0.22-µm are required to sterilize heat-sensitive media/solutions and to capture and culture viable test bacteria in control samples and eluates.
- 6.8 *Positive Displacement Pipette*—A pipette and pipette tips that accurately can dispense 10-µL volumes.
- 6.9 Refrigerator—A refrigerator at  $4 \pm 2^{\circ}$ C for storage of prepared culture media and reagents.
- 6.10 *Sterilizer*—Any suitable steam sterilizer capable of producing the conditions of sterilization is acceptable.
- 6.11 Timer(Stop-clock)—One that can be read for minutes and seconds.
- 6.11.1 Tap Water Temperature Regulator and Temperature Monitor—to monitor and regulate water temperature at 40  $\pm$  2°C.
- 6.11.2 *Water Faucet(s)*—to be located above the sink at a height that permits the hands to be held higher than the elbow during the washing procedure. Faucets with electronic sensors or those that are wrist-, elbow-, knee-, or foot-operated are preferred to avoid recontamination of the washed hands.

# 7. Materials and Reagents

- 7.1 Serological Pipettes—Sterile reusable or single-use pipettes of 10.0, 5.0, and 1.0-mL capacity.
- 7.2 Culture Plates—Petri plates of 100 mm diameter for culturing bacteria.

 ${\it Note}\ 1$ —Plastic culture ware may be purchased from most laboratory supply houses.

- 7.3 Culture Media and Supplements—Culture media and the types and ratios of supplements will vary depending on the type of test bacterium being used.
  - 7.4 Soil Load:
- 7.4.1 Fetal Bovine Serum, at a final concentration of 5 % in the bacterial inoculum.
  - 7.4.2 Tripartite Soil Load, as an alternative to serum.
- 7.4.2.1 Add 0.5 g of tryptone or yeast extract to 10 mL of phosphate buffer.
- 7.4.2.2 Add 0.5 g of bovine serum albumin (BSA) to 10 mL of phosphate buffer.
- 7.4.2.3 Add 0.04 g of bovine mucin to 10 mL of phosphate buffer.
- 7.4.3 Prepare the stock solutions separately and sterilize by passage through a 0.22  $\mu$ m pore diameter membrane filter, aliquot and store at either 4  $\pm$  2°C or –20  $\pm$  2°C; the stock solution of bovine mucin can be autoclave-sterilized.
- 7.4.4 To obtain a 500- $\mu$ L inoculum of the test inoculum, add to 340  $\mu$ L of the bacterial suspension, 35  $\mu$ L of tryptone or yeast extract (7.4.2.1), 25  $\mu$ L BSA (7.4.2.2), and 100  $\mu$ L mucin

- (7.4.2.3) stock solutions. This mixture contains approximately 2 g of total protein/L, which is roughly equivalent to the protein content of a 5 % solution of fetal bovine serum.
- 7.5 Standard Hard Water—The quality and disinfectant (for example, chlorine) residual in tap water can vary from site to site and also at different times at the same site. The use of standard hard water, therefore, is recommended here to avoid variations in results due to differences in tap water quality. Water prepared in accordance with AOAC 960.09 E and F (8) to a standard hardness of at least 200 ppm as calcium carbonate is used for dilution of the test substance, as the control solution to determine the baseline level of bacterial elimination, and to rinse the fingerpads after exposure to the test product. The standard hard water and tap water (if used) must first be tested to ensure that they do not have any activity against the test bacterium. If water with a different level of hardness is used for making the use-dilution of the test formulation, this change must be clearly justified and documented in the report.
- 7.6 *Test Agents*—At least two samples of the product shall be tested.
- 7.7 Diluent for Bacterial Titration—Normal saline (0.85 % NaCl) at pH 7.2 7.4, or appropriate buffer.
- 7.8 Eluent for Bacterial Recovery from Fingerpads—Normal saline or another suitable eluent.
- 7.9 *Plastic Vials*—Sterile screw-capped 2.0-mL vials with an inside diameter of about 8 mm will be required for demarcation of the fingerpads and to hold various test solutions.
- 7.10 *Miscellaneous Laboratory Ware*—Automatic pipettes, pipette tips, plastic vials for storing stock cultures.
- 7.11 *Broth*—Tryptose phosphate broth (TPB) or equivalent, to prepare the inoculum of the test organisms.
- 7.12 Agar—Trypticase soy agar (TSA) or equivalent, to recover and count the colonies of the test organism in control and test samples. The addition of any neutralizer(s) in such recovery media must first be properly validated. The use of selective-differential media for the detection of the test bacteria in such studies is not recommended because cells stressed or injured after germicide exposure may not form colonies on such agars.

## 8. Test Bacteria

8.1 Appendix X1 contains a list of suggested test bacteria. (Warning—The application of microorganisms to the skin may involve a health risk. Prior to applying the test organism to the skin, the antibiotic sensitivity profile of the strain should be determined. If the applied organism causes an infection, the antibiotic sensitivity profile should be made available to the attending clinician.)

# 9. Preparation of Test Bacterial Suspension

9.1 The stock culture should be at least two, but no more than five, 24-h broth transfers from the original ATCC stock over a period of one week. To prepare the test bacterial



suspension, add 0.1 mL of a 24-h broth culture to 100 mL of TSB (7.11) and incubate for 24  $\pm$  4 h at the appropriate temperature.

- 9.2 Add soil load (see 7.4), if required.
- 9.3 Swirl, vortex or shake the test bacterial cell suspension before withdrawal of each aliquot. Assay the in-put suspension for CFU at the beginning and end of the use period. Do not use a test suspension for more than 8 hours on any given day. The titer of the test suspension as CFU may not vary more than  $\pm$  0.5  $\log_{10}$  over an 8-h period.

# 10. Subjects

- 10.1 Recruit a sufficient number of healthy adult human subjects who have no clinical evidence of dermatoses, open wounds, hangnails or other skin disorders. The number of subjects required for a trial is dependent on the number of treatments within a study.
- 10.2 It is the responsibility of the user of this test method to arrange the necessary clearance for the use of adult subjects for testing and to obtain informed and written consent from those selected for the study before starting the tests.
- 10.3 Instruct subjects to avoid contact with antimicrobial products, other than the test substance, for the duration of their involvement in the study and for at least one week prior to the first test. This restriction includes antimicrobial-containing antiperspirants, deodorants, shampoos, lotions, and soaps, also such materials as acids, bases and solvents. Bathing in biocide treated pools, hot tubs, spas should also be avoided. Subjects are to be provided with a kit of non-antimicrobial personal care products for exclusive use for the duration of the study and rubber gloves to be worn when contact with antimicrobials or harsh chemicals cannot be avoided.

# 11. Procedure

- 11.1 Fig. 1 shows the main steps in the performance of this test method.
- 11.1.1 The subject will wash his/her hands with a nonmedicated soap for at least 10 s, rinse, and then dry them thoroughly with a clean paper or cloth towel. This procedure reduces variability in the test results by removing accumulated oil and dirt from the hands. Place about 3 to 5 mL of 70 % ( $\frac{1}{2}$ ) ethanol in the palm of one of the washed hands and instruct the subject to rub it well over the entire surface of both hands until the alcohol and water have evaporated completely (*Step 1*).
- 11.1.2 Press a thumbpad or fingerpad over the mouth of an empty plastic vial (see 7.9) to demarcate the area to receive the test inoculum (*Step 2*).
- Note 2—The demarcation also serves as a guide in the placement of the plastic vial on the skin in the subsequent elution of the inoculum.
- 11.1.3 Using a positive displacement pipette, deposit 10  $\mu$ L of the bacterial suspension, with or without a soil load, at the center of each demarcated area (*Step 3*).
- 11.1.4 It is recommended that thumbpads be used to determine the numbers of viable bacteria placed in each demarcated area (Input Control). Once a thumbpad has been contaminated, do not allow the inoculum to dry, but immediately elute it in accordance with 11.1.9.

- Note 3—The eluted fingerpads must be decontaminated (11.1.10) immediately so that they could be used later in holding the vials for eluting the fingerpads, as shown in (*Step 5*).
- 11.1.5 Allow the inoculum on all fingerpads to become visibly dry (~20 to 30 min) under ambient conditions (*Step 4*).
- 11.1.6 To determine the number of bacteria remaining viable after this drying period, elute the bacteria from two randomly selected fingerpads in accordance with 11.1.9 (*Step 8*).
- 11.1.7 Expose the dried inoculum on the required number of randomly selected fingerpads, and control or reference solution by placing 1.0 mL of the in-use dilution of the test substance, control, or reference solution in a plastic vial (7.9). Place a contaminated fingerpad over the mouth of the vial and invert it. Allow the contents of the vial to remain in contact with the contaminated area for 10 to 30 s (Step 5) while subjecting the vial to 10 full inversions. For viscous formulations, invert the vial and keep its contents in contact with the contaminated area for 10 to 20 s without any inversions. When testing a handwash agent, scrape the fingerpad on the inside rim of the vial to recover as much of the fluid as possible (Step 6); such scraping simulates the friction often applied during handwashing. For handrubs, do not scrape the fingerpad; the fluid remaining on the skin should be allowed to remain in place and recovered during the elution process.
- 11.1.8 If the test method needs to simulate post-treatment rinsing of hands, expose two fingerpads to 1 mL of hard water for 5 to 10 s (*Step 7*). Eluate the bacteria from the fingerpads (*Step 8*).
- 11.1.9 For bacterial elution, place the contaminated area of the thumb/finger over the mouth of a plastic vial (see 7.9) containing 1 mL of the eluent/neutralizer. Invert the vial with the pad still over it, and allow the eluent to remain in contact with the inoculated contaminated area for 5 to 10 s. Invert the vial 20 times with the pad still in place. Repeat the soak and inversion step once more. Finally, turn the vial upright and scrape the pad against the inside rim of the vial to recover as much of the fluid as possible ( $Step\ 8$ ).
- 11.1.10 Decontaminate the fingerpads by pressing them for 2 to 3 min directly over tissue paper or paper towel soaked in 70 % ( $\frac{1}{2}$ ) ethanol ( $\frac{1}{2}$ ).
- 11.1.11 Instruct the subjects to then wash their hands thoroughly with soap and water and dry them well before leaving the test area.
- 11.1.12 Titrate the eluates and controls for viable bacteria using a minimum of three culture plates for each dilution tested. If titrations cannot be carried out within 3 to 4 h of collection, store samples overnight at 4 to 10°C.
- 11.2 Neutralization of any antibacterial activity in the eluates must be properly validated. Add the test bacteria to the neutralized eluate and titrate it along with appropriate controls to demonstrate that there is no detectable loss in bacterial viability after exposure to the neutralized eluates.

# 12. Enumeration of Bacteria in Control and Test Samples

12.1 Process the control and test samples using either spread plating or membrane filtration or spiroplating to detect CFU of

Step #	
(1) The panelist washes hands with non-germicidal soap and tap water and dries them with paper towel. About 3-5 mL of 70% (v/v) ethanol is placed on the hands and they are rubbed together till dry.	70% PHANOL
(2) Each digit is pressed against the mouth an empty plastic vial (8 mm inside diam.) to demarcate the target area.	
(3) $10~\mu L$ of the test bacterial suspension, with or without soil load, is placed at center of each demarcated area. Inoculum from the two thumbpads is eluted immediately (Step 8 below) to act as 'input' control.	
(4) Inoculum on fingerpads allowed to become visibly dry (20-25 minutes). Two randomly selected fingerpads are eluted immediately (Step 8) at the end of drying ('baseline' control).	
(5) Dried inoculum on at least two randomly selected fingerpads is exposed to 1 mL of test product or control fluid in a plastic vial for desired contact time, with specified number of full inversions.	
(6) Skin scraped against inside lip of vial to collect as much fluid as possible. For waterless handwash agents or to determine bacterial elimination after exposure to the product alone, fingerpads can be eluted (Step 8) without further treatment.	THE RESERVE TO THE PARTY OF THE
(7) To simulate post-treatment rinsing of hands, fingerpads are exposed to 1 mL of hard water for 5-10 seconds. The test bacteria can be eluted (Step 8) at this stage or after drying of hands. To determine bacterial removal after the drying of washed hands, they can be dried in air or with paper or cloth towel for specified time and the bacteria recovered from them.	
(8) To elute the bacteria, the digit is placed on the mouth of a plastic vial with 1 mL of eluent and subjected to 20 full inversions; skin is scraped against inside lip of vial to collect as much fluid as possible. The eluates and controls are titrated for viable bacteria and log <sub>10</sub> reductions calculated.	
(9) The panelist decontaminates the hands by pressing the inoculated areas for 2-3 minutes against a tissue or paper towel soaked in 70% ethanol. The panelist then washes hands thoroughly with soap and water and dries them well before leaving the test area.	

Note 1—In Step 6, scraping of fingerpad surface is not required when testing handrubs.

FIG. 1 Procedure For In Vivo Evaluation Of The Bacteria-Eliminating Activity Of Handwash And Handrub Agents
Using The Fingerpad Test

the test bacterium. The membrane filtration method (a) has the advantage of allowing larger sample volumes to be processed especially when the viable counts in the eluates may be low, and (b) permits the rinsing of the filters to reduce microbicidal residues, which may interfere with bacterial recovery. The pour-plate technique is not recommended for S. marcescens because subsurface colonies may not exhibit the characteristic red pigment.

12.2 Incubate the inoculated plates and sterility controls of the recovery medium for 48  $\pm$  4 h at the appropriate temperature.

12.3 Count colonies and convert data to  $log_{10}$  units.

# 13. Repetitions and Statistical Evaluations

- 13.1 Preform the testing at least twice with a minimum of three subjects for each formulation tested.
- 13.2 This test method is designed to include the two thumbpads to determine the viable bacteria placed on the fingerpads (input control), two fingerpads to assess the number of bacteria remaining viable after the drying of the inoculum (dry control), two fingerpads to determine the extent of bacterial elimination after treatment with standard hard water



alone, and four fingerpads to assess the level of bacterial elimination after the combined effect of exposure to the test product and the post-treatment rinse with the standard hard water. For waterless agents, the hard water rinse is not called for.

13.3 The difference in the number of viable bacteria in the inoculum control and the dry control represents the loss in bacterial viability due to the drying of the inoculum. The number of viable bacteria remaining after the drying of the inoculum must be used as the baseline to determine the extent of bacterial elimination after treatment with the control solution or the test product.

### 14. Precision and Bias

14.1 *Precision*—The precision of this test method within two laboratories has been determined and compared with that

of a whole-hand protocol (6). The efficiency of the test method of bacterial elution from the fingerpads has been tested with two different types of bacteria and has been found to be about 90 %.

# 15. Keywords

15.1 bacteria-eliminating activity; eluent; *Escherichia coli*; fingerpads; germicidal soap; hygienic handwashing; infection control; *in vivo* testing; *Serratia marcescens*; skin flora; soil load; standard hard water; *Staphylococcus aureus*; *Staphylococcus epidermidis* 

### **APPENDIX**

(Nonmandatory Information)

### X1. BACTERIA RECOMMENDED FOR USE IN THIS TEST

- X1.1 The selection of the following test bacteria is based on their (a) relative safety to the subjects as well as experimenters, (b) ability to grow to titers sufficiently high for testing, (c) ease of recovery and identification in the laboratory, and (d) potential to spread through contaminated hands. Other strains or types of bacteria may be substituted provided they meet the preceding criteria.
- X1.1.1 Serratia marcescens (ATCC 14756). A Gramnegative bacillus having stable pigmentation at 25°C.
- X1.1.2 Escherichia coli (ATCC 11229), A Gram-negative bacillus often used as a representative enteric bacterium.
- X1.1.3 Acinetobacter baumannii (ATCC 19606), A Gramnegative normal inhabitant of human skin, but a significant opportunistic pathogen with relatively high resistance to many antibiotics and microbicides.
- X1.1.4 Staphylococcus aureus (ATCC 6538). A normal inhabitant of human body and also an important nosocomial pathogen.
- X1.1.5 Staphylococcus epidermidis (ATCC 14990), A Gram-positive test organism that is considered safe for the experimental contamination of the skin of adults.

#### REFERENCES

- Draft CDC Guideline for Hand Hygiene in Healthcare Settings, U.S. Department of Health and Human Services, U.S. Government Printing Office, Washington, DC, 2001.
- (2) Ansari, S. A., Sattar, S. A., Springthorpe, V. S., Wells, G. A., and Tostowaryk, W., "Rotavirus Survival on Human Hands and Transfer of Infectious Virus Animate and Nonporous Inanimate Surfaces," *Journal of Clinical Microbiology*, Vol 26, 1988, pp. 1513–1518.
- (3) Mbithi, J. N., Springthorpe, V. S., Boulet, J., and Sattar, S. A., "Survival of Hepatitis A Virus on Human Hands and Its Transfer on Contact with Animate and Inanimate Surfaces," *Journal of Clinical Microbiology*, Vol 30, 1992, pp. 757–763.
- (4) Sattar, S. A., Abebe, M., Bueti, A., Jampani, H., and Newman, J. "Determination of the activity of an alcohol-based hand gel against human adeno-, rhino-, and rotaviruses using the fingerpad method," *Infection Control & Hospital Epidemiology*, Vol 21, 2000, pp. 516–519.

- (5) CDC-NIH, Biosafety in Microbiological and Biomedical Laboratories, 4th ed., U.S. Department of Health and Human Services, U.S. Government Printing Office, Washington, DC, 1999.
- (6) Ansari, S. A., Sattar, S. A., Springthorpe, V. S., Wells, G. A., and Tostowaryk, W., "In Vivo Protocol for Testing Efficacy of Hand-Washing Agents Against Viruses and Bacteria: Experiments With Rotavirus and Escherichia coli," Applied and Environmental Microbiology, Vol 55, 1989, pp. 3113–3118.
- (7) Ansari, S. A., Springthorpe, V. S., and Sattar, S. A., "Comparison of Cloth-, Paper- and Warm Air-Drying in Eliminating Viruses and Bacteria from Washed Hands," *American Journal of Infection Control*, Vol 19, 1991, pp. 243–249.
- (8) AOAC International, Official Methods of Analysis of the AOAC, Arlington, VA, 1990.



ASTM International takes no position respecting the validity of any patent rights asserted in connection with any item mentioned in this standard. Users of this standard are expressly advised that determination of the validity of any such patent rights, and the risk of infringement of such rights, are entirely their own responsibility.

This standard is subject to revision at any time by the responsible technical committee and must be reviewed every five years and if not revised, either reapproved or withdrawn. Your comments are invited either for revision of this standard or for additional standards and should be addressed to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee, which you may attend. If you feel that your comments have not received a fair hearing you should make your views known to the ASTM Committee on Standards, at the address shown below.

This standard is copyrighted by ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, PA 19428-2959, United States. Individual reprints (single or multiple copies) of this standard may be obtained by contacting ASTM at the above address or at 610-832-9585 (phone), 610-832-9555 (fax), or service@astm.org (e-mail); or through the ASTM website (www.astm.org). Permission rights to photocopy the standard may also be secured from the Copyright Clearance Center, 222 Rosewood Drive, Danvers, MA 01923, Tel: (978) 646-2600; http://www.copyright.com/