INTERNATIONAL STANDARD

ISO 8578

First edition 1997-04-15

Optics and optical instruments — Microscopes — Marking of objectives and eyepieces

Optique et instruments d'optique — Microscopes — Marquage des objectifs et des oculaires



ISO 8578:1997(E)

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

International Standard ISO 8578 was prepared by Technical Committee ISO/TC 172, Optics and optical instruments, Subcommittee SC 5, Microscopes and endoscopes.

Annex A forms an integral part of this International Standard.

© ISO 1997

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from the publisher.

International Organization for Standardization
Case postale 56 • CH-1211 Genève 20 • Switzerland
Internet central @iso.ch
X.400 c=ch; a=400net; p=iso; o=isocs; s=central

Printed in Switzerland

Introduction

This International Standard has been prepared in order to define clearly the data relating to optical characteristics with which microscope objectives and eyepieces have to be marked and the positioning of such marking on the component to enable correct use of the microscope. In addition to data which have to be marked, recommendations for the marking of additional information relating to several other optical characteristics are given.



Optics and optical instruments — Microscopes — Marking of objectives and eyepieces

1 Scope

This International Standard specifies the format for the marking of data for optical characteristics on microscope objectives and eyepieces and the positioning of this data. It makes recommendations for the marking of additional information, particularly colour coding of rings designating the magnification of objectives and the immersion media with which they are used.

2 Objectives

2.1 Mandatory markings on objectives

The markings on objectives shall be as given in table 1.

2.2 Recommended additional markings on objectives

The marking of additional data is optional. If additional markings are used, they should be as given in table 2.

2.3 Recommendation for arrangement of the marking

It is recommended that the markings in column A of table 3 should be placed above or before the markings of column B, which in turn should be placed above or before those of column C.

3 Eyepieces

3.1 Mandatory markings on eyepieces

The markings on eyepieces shall be as specified in table 4.

3.2 Recommended additional markings on eyepieces

The marking of additional data is optional. If additional markings are used, they should be as given in table 5.

ISO 8578:1997(E) © ISO

Table 1 — Mandatory markings on objectives

Optical property	Feature to be marked	Example of marking ¹⁾	Remarks	
Magnification	Lateral magnification of objectives for a finite image distance	100	Magnification and numerical aperture should be separated by an oblique stroke, e.g. 100/1,30	
	Lateral magnification of objectives for an infinite image distance	100×	The marked value of the magnification of infinity-corrected objectives is valid only in combination with the related tube lens. The marking of the symbol "x" has been introduced as an additional designation of magnification on infinity-corrected objectives	
Aperture	Numerical aperture	/1,30	The numerical aperture shall be stated to at least 2 decimal places	
Immersion medium	Oil for immersion oil	OIL	For additional marking, colour-coded rings can be used (see 2.2)	
	W for water	W		
	Glyc for glycerol	GLYC		
	Other		The requirement to use any other immersion medium shall be indicated	
Tube length	Length, in mm, for objective of finite primary image distance	160	Tube length and cover glass thickness shall be separated by an oblique stroke, e.g. 160/0,17, 160/–, 160/0, and shall be printed in a smaller size than the data for magnification and aperture	
	Symbol ∞ for objective of primary image distance ∞	∞		
Equivalent cover glass thickness	Thickness, in mm	/0	For objectives that are corrected for use with uncovered specimens only, the figure "0" shall be indicated after the oblique stroke	
		/0,17	For objectives that are corrected for use only with a cover glass, the value of the cover glass thickness to be used shall be indicated, in millimetres, after the oblique stroke, e.g. 0,17	
		/-	For objectives that can be used without a cover glass or with a cover glass up to 0,17 mm thickness, the symbol "—" shall be positioned after the oblique stroke	
Phase contrast	Symbol PH	PH 2	A figure after the symbol indicates the associated annular diaphragm	
Polarizing microscopy system	Symbol POL	POL		
Flatness of field	Symbol PLAN or PL	PLAN	The symbols "PLAN" or "PL" shall only be used if there is a minimum flat field of 18 mm diameter when the objective is used with a tube factor of 1×	

Table 1 (concluded)

Feature to be marked	Example of marking ¹⁾	Remarks
Achromat		Achromatic objectives require no marking to indicate the nature of their chromatic correction
Apochromat, symbol APO	APO	Objectives with a chromatic correction intermediate between achromat and apochromat are generally marked with the manufacturer's designation indicating such correction
Limiting values of numerical aperture	/1,30–0,8	The lower and upper limits of the numerical aperture range controlled by the iris diaphragm shall be marked at the position where the value for the numerical aperture is usually marked
Name or symbol of identification		
	Achromat Apochromat, symbol APO Limiting values of numerical aperture Name or symbol of	Achromat Apochromat, symbol APO Limiting values of numerical aperture Apochromat, symbol APO APO /1,30–0,8

ISO 8578:1997(E) © ISO

Table 2 — Recommended additional markings on objectives

Optical property	Feature to	be marked	Example of marking ¹⁾	Remarks
Magnification	Value	Colour of ring		
	1/1,25	Black		
	1,6/2	Grey		
	2,5/3,2	Brown		
	4/5	Red		
	6,3/8	Orange		
	10/12,5	Yellow		
	16/20	Light green		
	25/32	Dark green		
	40/50	Light blue		
	63/80	Dark blue		
	100 125 160	White		
Immersion medium	Medium	Colour of ring		To avoid confusion, it is recommended that a coloured ring indicating the immersion medium should only be used in conjunction with a second ring indicating the magnification
	Air	None		
	Oil	Black		
	Water	White		
	Glyc	Orange		
	Others	Red		
Phase contrast	The entire apart from coloured rin manufactur	the ngs and the er's name,		The marking of the manufacturer's name may be in any colour
Polarizing microscopy systems	The entire is apart from to coloured ring manufactures shall be in a	the ngs and the er's name,		The marking of the manufacturer's name may be in any colour
Differential interference contrast	Symbol DIC		DIC	
Objectives for Epi illumination	Symbol EP	I	EPI	

Table 2 (concluded)

Optical property	Feature to be marked	Example of marking ¹⁾	Remarks
Objectives for Epi illumination, brightfield and darkfield	Symbol D	D	The symbol EPI may be marked in addition
Long working distance	Symbol L	L	
Country of manufacture			The marking of the country of origin is mandatory in several countries

Table 3 — Recommendation of arrangement of markings on objectives

A 1)	B 1)	C 1)
Flatness of field	Magnification	Immersion medium
State of chromatic correction	Numerical aperture	Phase contrast
Long working distance		Polarizing microscopy system
		Differential interference contrast
		Objectives for brightfield and darkfield Epi illumination

¹⁾ A to precede B, to precede C.

NOTES

- 1 If an additional coloured ring is used in accordance with table 2 to identify the immersion medium, this ring should be placed closer to the front lens than the coloured ring used to indicate the magnification.
- 2 The marking relating to tube length and cover thickness, specified in table 1, can be in the position of the markings given in either column A or column C.

Table 4 — Mandatory markings on eyepieces

Optical property	Feature to be marked	Example of marking	Remarks
Magnification	Visual magnification	10×	Visual magnification and field-of-view number shall be separated by an oblique stroke, e.g. 10×/18
Field of view	Diameter, in mm	/18	
Manufacturer	Name or symbol of identification		

ISO 8578:1997(E) © ISO

Table 5 — Recommended additional markings on eyepieces

Optical property	Feature to be marked	Example of marking ¹⁾	Remarks
Suitable for spectacle wearers	Symbol &	<i>6</i> 67	Marking in conjunction with the magnification and field-of-view number, for example 10×/18 &
Centred and aligned crosslines	Symbol ⊕	⊕	Only for eyepieces with centred crosslines aligned with an orienting pin
Type of correction	Compensating symbol C	С	If the eyepiece is corrected to project a real image, the additional marking "PHOT" shall
	Flatfield symbol PL	PL	follow the indication of the type of correction
1) Capital or lower case	e letters optional.		

Annex A

(normative)

Validity of other markings

A.1 General

The purpose of this International Standard is to specify one type of marking or symbol to clearly define each optical characteristic of objectives and eyepieces.

However, in order to give the manufacturers enough time to adapt their entire production to this International Standard, a transition period is necessary.

"Phaco"

For a transition period of 10 years, the symbols in A.2 and A.3 shall also be allowed.

NOTE — Symbols may be given in lower case or capital letters.

A.2 Mandatory data on objectives (see table 1)

Phase contrast:

Polarizing microscopy system: "P" or "PO"

A.3 Recommended additional data on objectives (see table 2)

For objectives for EPI illumination: "M"

For objectives for EPI illumination, brightfield and "BD" or "HD"

darkfield:

For objectives with long working distance: "LD" or "LWD"

ICS 37.020

Descriptors: optics, optical equipment, microscopes, optical microscopes, objectives, eyepieces, marking.

Price based on 7 pages