TECHNICAL SPECIFICATION

150/TS 11869

IDF/RM 150

First edition 2012-01-15

Fermented milks — Determination of titratable acidity — Potentiometric method

Laits fermentés — Détermination de l'acidité titrable — Méthode potentiométrique





ISO/TS 11869:2012(E) IDF/RM 150:2012(E)



COPYRIGHT PROTECTED DOCUMENT

© ISO and IDF 2012

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO or IDF at the respective address below.

ISO copyright office Case postale 56 • CH-1211 Geneva 20 Tel. + 41 22 749 01 11 Fax + 41 22 749 09 47 E-mail copyright@iso.org Web www.iso.org

Published in Switzerland

International Dairy Federation Silver Building • Boulevard Auguste Reyers 70/B • B-1030 Brussels Tel. + 32 2 733 98 88 Fax + 32 2 733 04 13

E-mail info@fil-idf.org Web www.fil-idf.org

Foreword

The International Organization for Standardization (ISO) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

In other circumstances, particularly when there is an urgent market requirement for such documents, a technical committee may decide to publish other types of document:

- an ISO Publicly Available Specification (ISO/PAS) represents an agreement between technical experts in an ISO working group and is accepted for publication if it is approved by more than 50 % of the members of the parent committee casting a vote;
- an ISO Technical Specification (ISO/TS) represents an agreement between the members of a technical committee and is accepted for publication if it is approved by 2/3 of the members of the committee casting a vote.

An ISO/PAS or ISO/TS is reviewed after three years in order to decide whether it will be confirmed for a further three years, revised to become an International Standard, or withdrawn. If the ISO/PAS or ISO/TS is confirmed, it is reviewed again after a further three years, at which time it must either be transformed into an International Standard or be withdrawn.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO/TS 11869|IDF/RM 150 was prepared by Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 5, *Milk and milk products*, and the International Dairy Federation (IDF). It is being published jointly by ISO and IDF.

This first edition cancels and replaces ISO 11869:1997, which has been technically revised.

ISO/TS 11869:2012(E) IDF/RM 150:2012(E)

Foreword

IDF (the International Dairy Federation) is a non-profit organization representing the dairy sector worldwide. IDF membership comprises National Committees in every member country as well as regional dairy associations having signed a formal agreement on cooperation with IDF. All members of IDF have the right to be represented on the IDF Standing Committees carrying out the technical work. IDF collaborates with ISO in the development of standard methods of analysis and sampling for milk and milk products.

The main task of Standing Committees is to prepare International Standards. Draft International Standards adopted by the Standing Committees are circulated to the National Committees for endorsement prior to publication as an International Standard. Publication as an International Standard requires approval by at least 50 % of IDF National Committees casting a vote.

In other circumstances, particularly when there is an urgent market requirement for such documents, a Standing Committee may decide to publish an other type of normative document which is called by IDF: Reviewed method. Such a method represents an agreement between the members of a Standing Committee and is accepted for publication if it is approved by at least 50 % of each National Committee representative present casting a vote. A Reviewed method is equal to an ISO/PAS or ISO/TS and will, therefore, also be published jointly under ISO conditions.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. IDF shall not be held responsible for identifying any or all such patent rights.

ISO/TS 11869|IDF/RM 150 was prepared by the International Dairy Federation (IDF) and Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 5, *Milk and milk products* and the International Dairy Federation (IDF). It is being published jointly by ISO and IDF.

All work was carried out by the Joint ISO-IDF Project Group on *Titratable acidity* of the Standing Committee on *Analytical methods for composition* under the aegis of its project leader, Dr. J. Floor (ZA).

This edition of ISO/TS 11869|IDF/RM 150 cancels and replaces IDF 150:1991, of which it constitutes a technical revision.

Fermented milks — Determination of titratable acidity — Potentiometric method

1 Scope

This Technical Specification specifies a potentiometric method for the determination of the titratable acidity of natural yoghurt, flavoured yoghurt, fruit yoghurt, drinking yoghurt, fresh cheese with or without fruit, buttermilk with or without fruit, and other fermented milk products.

2 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

2.1

titratable acidity of fermented milks

amount in millilitres of a 0,1 mol/l sodium hydroxide solution required to titrate 10 g of product to pH $8,30\pm0,01$

NOTE The titratable acidity is expressed in millimoles per 100 g.

3 Principle

A test portion is suspended in water. The suspension is titrated potentiometrically against sodium hydroxide solution [c(NaOH) = 0.1 mol/l] to pH 8,30 \pm 0,01. The titratable acidity is calculated.

4 Reagents

Use only reagents of recognized analytical grade, unless otherwise specified, and only distilled or deionized water, freed from carbon dioxide by boiling for 10 min before use.

4.1 Sodium hydroxide, standard volumetric solution, $c(NaOH) = 0.1 \text{ mol/l} \pm 0.002 \text{ mol/l}$, carbonate free.

Protect this solution against absorption of carbon dioxide (CO₂) either by connecting a washing bottle with 10 % sodium hydroxide solution to the burette which itself contains the sodium hydroxide solution or by connecting a small tube with fresh sodium hydroxide or calcium oxide to the end of the burette to obtain a closed system.

NOTE CO₂ is either bound in the washing bottle or in the tube to protect the solution in the burette against absorption which would influence the concentration.

ISO/TS 11869:2012(E) IDF/RM 150:2012(E)

5 Apparatus

Usual laboratory equipment and, in particular, the following.

- **5.1 Analytical balance**, capable of being read to the nearest 0,01 g.
- **5.2 pH meter**, correctly calibrated in the range pH 7 to 10 following normal laboratory pH-calibration procedures.
- 5.3 Spoon or spatula.
- **5.4** Homogenizer, e.g. a macerator [Ultra-Turrax¹⁾ or equivalent].
- **5.5 Burette**, of capacity 25 ml or 50 ml, graduated at least at every 0,05 ml, ISO 385^[1], class A. Alternatively, an **automatic burette**, ISO 8655-3^[4], fulfilling the same requirements may be used.

NOTE Instead of carrying out a manual titration, it is also possible to use an automatic titrator.

5.6 Water bath, capable of maintaining a temperature of 38 °C \pm 1 °C.

6 Sampling

Sampling is not part of the method specified in this Technical Specification. A recommended sampling method is given in ISO $707[IDF\ 50^{[2]}]$.

It is important the laboratory receive a truly representative sample which has not been damaged or changed during transport or storage.

7 Preparation of test sample

7.1 Natural yoghurt, flavoured yoghurt, drinking yoghurt and other fermented milks

Bring the sample to a temperature of 22 °C \pm 2 °C. Mix the sample carefully by means of a spoon or spatula (5.3) or homogenizer (5.4), using a rotary motion which passes from the lower layers to the surface layers of the sample so as to displace and mix them well.

7.2 Fruit yoghurt and other fermented milk products with added fruit

Bring the sample to a temperature of 22 °C \pm 2 °C. Homogenize it using an appropriate device (5.4), in order to facilitate the grinding and dispersion of the fruit.

If fat separation is observed in the sample, raise the temperature of the sample in the water bath (5.6) to $38 \,^{\circ}\text{C} \pm 1 \,^{\circ}\text{C}$ for a better homogenization. Thereafter, cool the sample again to a temperature of 22 $\,^{\circ}\text{C} \pm 2 \,^{\circ}\text{C}$.

¹⁾ Ultra-Turrax is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO or IDF of this product.

8 Procedure

8.1 Test portion

Weigh, in a 50 ml beaker, approximately 10 g of prepared test sample (Clause 7), to the nearest 0,01 g. Add approximately 10 ml of water and mix.

8.2 Determination

- **8.2.1** Insert the electrode of the pH meter (5.2) into the suspension (8.1).
- **8.2.2** Use the sodium hydroxide solution (4.1), while stirring, to titrate the contents of the beaker until the pH remains stable at $8,30 \pm 0,01$ for 4 s to 5 s.

If an automatic titrator is used, it shall provide for such a stop delay time as well.

Record the volume, in millilitres, of sodium hydroxide solution used, at least to the nearest 0,05 ml.

9 Calculation and expression of results

9.1 Calculation

Calculate the titratable acidity, *I*, in millimoles of sodium hydroxide per 100 g, using the following equation:

$$I = \frac{V \times 10}{m}$$

where

- V is the volume, in millilitres, of the sodium hydroxide solution (4.1) used for the titration (8.2.2);
- m is the mass, in grams, of the test portion (8.1).

NOTE For the expression of results in grams of lactic acid per 100 g of product, see Annex A.

9.2 Expression of test results

Express the test results to two decimal places.

10 Precision

10.1 Interlaboratory test

The method has been tested on various fermented milk products in six different laboratories. Each laboratory obtained the products locally. See Reference [5]. A summary of the results is given in Annex B. The value for repeatability derived from this interlaboratory study was calculated in accordance with ISO 5725-1^[3].

10.2 Repeatability

The absolute difference between two independent single test results, obtained using the same method on identical test material in the same laboratory by the same operator using the same equipment within a short interval of time, will in not more than 5 % of cases be greater than 0,20 mmol/100 g.

ISO/TS 11869:2012(E) IDF/RM 150:2012(E)

Test report

The test report shall contain at least the following information:

- all information necessary for the complete identification of the sample; a)
- the sampling method used, if known; b)
- the method used, with reference to this Technical Specification (ISO/TS 11869|IDF/RM 150:2012); c)
- all operating details not specified in this Technical Specification, or regarded as optional, together with d) details of any incidents which may have influenced the test result(s);
- the test result(s) obtained; e)
- f) if the repeatability has been checked, the final quoted result obtained.

Annex A (informative)

Alternative expression of test results

In the previous version of this method, test results were expressed in grams of lactic acid per 100 g of product, w. This is slightly misleading as the method is not a procedure for the quantitative determination of lactic acid. However, if expression of the results in this way is preferred, the formula for the calculation of results is as follows:

$$w = \frac{V \times 0.9}{m}$$

where

V is the volume, in millilitres, of the sodium hydroxide solution (4.1) used for the titration (8.2.2);

m is the mass, in grams, of the test portion (8.1).

Annex B (informative)

Interlaboratory trial


An interlaboratory trial involving six laboratories was carried out^[5]. Each laboratory locally obtained five or six samples of various cultured products and analysed the samples in six-fold for titratable acidity using the method specified in this Technical Specification. The results are summarized in Table B.1.

Table B.1 — Results of interlaboratory trial for determination of titratable acidity of fermented milks

Natural yoghurt	Laboratory	Product	Average result mmol/100 g	Standard deviation mmol/100 g	
Flavoured yoghurt	1	Natural yoghurt	12,57	0,080	
Natural buttermilk 8,17 0.033		Fruit yoghurt	9,53	0,058	
Natural buttermilk 8,17 0.033		Flavoured yoghurt	11,28	0,027	
Flavoured drinking yoghurt 8,59 0,008		Natural buttermilk	8,17	0.033	
Natural yoghurt 8,62 0,058		Flavoured drinking yoghurt	10,03	0,021	
Natural (bio)yoghurt 12,85 0,034		Flavoured drinking yoghurt	8,59	0,008	
Buttermilk with fruit	2	Natural yoghurt	8,62	0,058	
Natural fermented milk 10,87 0,042 Natural fermented milk 9,48 0,038 Natural yoghurt 9,90 0,061 Fruit yoghurt 9,25 0,051 Flavoured yoghurt 10,10 0,038 Drinking yoghurt with fruit 10,32 0,027 Natural buttermilk 9,20 0,099 Buttermilk with fruit 9,92 0,012 Natural yoghurt 11,29 0,193 Flavoured drinking yoghurt 9,38 0,147 Flavoured drinking yoghurt 10,57 0,082 Flavoured (bio)yoghurt 10,85 0,152 Natural yoghurt 11,00 0,093 Natural yoghurt 11,00 0,093 Natural yoghurt 8,78 0,113 Fruit yoghurt 10,52 0,090 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Natural (bio)yoghurt	12,85	0,034	
Natural fermented milk 9,48 0,038		Buttermilk with fruit	8,85	0,180	
Natural yoghurt 9,90 0,061		Natural fermented milk	10,87	0,042	
Fruit yoghurt 9,25 0,051 Flavoured yoghurt 10,10 0,038 Drinking yoghurt with fruit 10,32 0,027 Natural buttermilk 9,20 0,099 Buttermilk with fruit 9,92 0,012 Natural yoghurt 11,29 0,193 Flavoured yoghurt 9,38 0,147 Flavoured drinking yoghurt 10,57 0,104 Natural fermented milk 10,85 0,152 Natural yoghurt 11,00 0,093 Natural yoghurt 11,00 0,093 Natural yoghurt 8,78 0,113 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 10,52 0,090 Flavoured yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Natural fermented milk	9,48	0,038	
Flavoured yoghurt	3	Natural yoghurt	9,90	0,061	
Drinking yoghurt with fruit 10,32 0,027 Natural buttermilk 9,20 0,099 Buttermilk with fruit 9,92 0,012 Natural yoghurt 11,29 0,193 Flavoured yoghurt 9,38 0,147 4		Fruit yoghurt	9,25	0,051	
Drinking yoghurt with fruit 10,32 0,027 Natural buttermilk 9,20 0,099 Buttermilk with fruit 9,92 0,012 Natural yoghurt 11,29 0,193 Flavoured yoghurt 9,38 0,147 4 Flavoured drinking yoghurt 8,79 0,082 Flavoured (bio)yoghurt 10,57 0,104 Natural fermented milk 10,85 0,152 Natural yoghurt 11,00 0,093 Natural yoghurt 11,00 0,093 Natural yoghurt 8,78 0,113 5 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Flavoured yoghurt	10,10	0,038	
Buttermilk with fruit 9,92 0,012 Natural yoghurt 11,29 0,193 Flavoured yoghurt 9,38 0,147 4		Drinking yoghurt with fruit	10,32	0,027	
Natural yoghurt		Natural buttermilk	9,20	0,099	
Flavoured yoghurt 9,38 0,147 Flavoured drinking yoghurt 8,79 0,082 Flavoured (bio)yoghurt 10,57 0,104 Natural fermented milk 10,85 0,152 Natural yoghurt 11,00 0,093 Natural yoghurt 8,78 0,113 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Buttermilk with fruit	9,92	0,012	
Flavoured drinking yoghurt 8,79 0,082	4	Natural yoghurt	11,29	0,193	
Flavoured (bio)yoghurt		Flavoured yoghurt	9,38	0,147	
Natural fermented milk 10,85 0,152 Natural yoghurt 11,00 0,093 Natural yoghurt high fat 8,78 0,113 5 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Flavoured drinking yoghurt	8,79	0,082	
Natural yoghurt 11,00 0,093 Natural yoghurt high fat 8,78 0,113 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Flavoured (bio)yoghurt	10,57	0,104	
Natural yoghurt high fat 8,78 0,113 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Natural fermented milk	10,85	0,152	
5 Fruit yoghurt 8,99 0,059 Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053	5	Natural yoghurt	11,00	0,093	
Fruit yoghurt 10,52 0,090 Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Natural yoghurt high fat	8,78	0,113	
Flavoured yoghurt 8,71 0,053 Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Fruit yoghurt	8,99	0,059	
Fruit yoghurt 10,91 0,046 Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Fruit yoghurt	10,52	0,090	
Natural fresh cheese 11,01 0,067 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Flavoured yoghurt	8,71	0,053	
6 Natural fresh cheese 10,46 0,063 Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053	6	Fruit yoghurt	10,91	0,046	
Fresh cheese with fruit 10,77 0,039 Fresh cheese with fruit 11,77 0,053		Natural fresh cheese	11,01	0,067	
Fresh cheese with fruit 11,77 0,053		Natural fresh cheese	10,46	0,063	
		Fresh cheese with fruit	10,77	0,039	
Grand mean repeatability standard deviation s_r 0,069		Fresh cheese with fruit	11,77	0,053	
	Grand mean re	Grand mean repeatability standard deviation s_r			
Repeatability limit $r = 2.8 s_r$ 0,20	Repeatability lin	Repeatability limit $r = 2.8 s_r$			

Bibliography

- [1] ISO 385, Laboratory glassware Burettes
- [2] ISO 707|IDF 50, Milk and milk products Guidance on sampling
- [3] ISO 5725-1, Accuracy (trueness and precision) of measurement methods and results Part 1: General principles and definitions
- [4] ISO 8655-3, Piston-operated volumetric apparatus Part 3: Piston burettes
- [5] Interlaboratory trial on the potentiometric method for the determination of titratable acidity in fermented milk products according to ISO/TS 11869|IDF/RM 150. *Bull. IDF* in press



ICS 67.100.99

Price based on 7 pages